SightLife Surgical )™

# **DONOR TREPHINE VACUUM PUNCH**

## **PRODUCT INFORMATION AND INSTRUCTIONS**

**SightLife Surgical**

101 North Chestnut Street, Ste 303  
Winston-Salem, NC 27101

[www.sightlifesurgical.com](http://www.sightlifesurgical.com)  
Tél. : +1 (844) 526 7632

EN

## SYMBOLS

**LOT** Lot Number

 Do not reuse

**STERILE EO**

**REF** Catalog Number

Sterilised by Ethylene Oxide

 Use Until Date

 Do not use if product is opened or damaged

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
## Instructions For Use

### Description.

The Donor Trephine Vacuum Punch is a single use only punch and trephine system for the trephination of the donor cornea button. The system consists of a vacuum base, a trephine guide and a standard SightLife Surgical trephine (see Fig. 1). It also contains a protector cover as an accessory.

### Intended Use.

The Donor Trephine Vacuum Punch is designed to allow the donor corneal tissue to be held in the most relaxed state in the bowl, while a vacuum is drawn. The vacuum holds the tissue gently in the optimal position, in preparation for the cutting of the donor button. The trephine, supplied as part of the system, is ultra sharp and makes a clean vertical cut through the corneal tissue.

### CAUTION

- This device is intended for SINGLE USE ONLY. It is not intended for repeat use.
- The vacuum punch is supplied sterile and ready for use. Do not use if the peel pouch has been opened or damaged.
- The **SightLife Surgical Trephine** has an ultra sharp cutting blade. Caution should be exercised at all times when handling and using the trephine.
- The **SightLife Surgical Trephine** is supplied with the blade in perfect condition. Should the device be accidentally dropped onto a hard surface prior to use, the blade integrity should be checked under the operating microscope. If found to be damaged, the trephine must be discarded.
- It is recommended that the trephine be left in the packaging blister until required, to ensure that the blade remains in perfect condition.

### Surgical Technique.

1. Remove the **punch base** and syringe from the packaging.
2. Remove the cornea from the storage medium by grasping the rim of the surrounding sclera with toothed forceps. Dry the anterior surface with a lint free PVA sponge spear.
3. In preparation for the procedure, it is suggested that an assistant depresses the plunger of the syringe prior to the placement of the donor tissue in the bowl of the vacuum base.
4. Position the corneal tissue in the bowl of the base.
5. Confirm the position of the corneal tissue and if necessary use a micro hook to manipulate and position the tissue optimally.

6. Release the plunger on the syringe to create the vacuum.
7. Once the vacuum has been established, position the **trephine guide** over the bowl so that it will drop vertically onto the base. Press down gently until the two parts are aligned correctly.
8. Remove the trephine from the packaging. Position the trephine on the guide and gently lower the blade onto the endothelial surface.
9. Press down on the trephine with gentle but firm pressure to push the blade through the corneal tissue until the blade comes to a complete stop. Rotate the trephine carefully through 10 degrees to ensure that the cut is complete.
10. Carefully remove the trephine and the trephine guide together as one piece. The button will remain in the bowl, held gently in place by the vacuum.
11. In the unlikely event that the corneal button remains in the lumen of the trephine, a balanced salt solution or a viscoelastic substance should be gently dripped into the centre of the blade to dislodge the button.
12. Slightly depress the plunger of the syringe or cut the silicone tube to break the vacuum.
13. Cover the corneal tissue with the protector cap.
14. When the recipient cornea has been prepared, the protector cap may be removed and the corneal button removed with a corneal spoon or micro forceps.

## **Sterilisation.**

### DO NOT RE-STERILISE.

The system is supplied sterile by Ethylene Oxide and is ready for use. SightLife Surgical assumes no liability for devices that have been re-sterilised by healthcare facilities.

## **Disposal.**

For safe disposal, the trephine should be considered as 'sharp' and must be immediately placed in a secure puncture-resistant bin (conforming to BS7320) which is suitable for incineration. All the components should be disposed of in accordance with hospital protocol.

## **HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:**

1. Single use devices have not been validated for re-use.  
If you re-use a device you may be held **Legally Liable** for the safe performance.

2. Cross-contamination and infection risks to patients . Including transmission of:
  - CJD & Variant CJD.
  - Prion Diseases.
  - Bacterial Endotoxins.
  - Hepatitis B & Hepatitis C.
  - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:
  - **Plastics:** Can be weakened, warped or become brittle.
  - **Metals:** Can be damaged or subjected to rusting. Blade edges will be damaged or blunt. Needles may be bent from the original shape.
  - **Adhesives:** Can degrade if re-sterilised causing parts to dislodge.
  - **Other materials:** May degrade, becoming unacceptable when compared to original manufacturing criteria.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

## PRODUCT DETAILS

SightLife Surgical Donor Trehphine Vacuum Punch  
Sizes 6.0mm-9.75mm in 0.25mm increments

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