EN

SYMBOLS

LOT

Lot Number

REF

Catalog Number

STERILE EO

Sterilised by Ethylene Oxide

Use Until Date

Do not use if product is opened or damaged

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

clockwise =
de

Marker
pin

SightLife Surgical RECIPIENT VACUUM TREPHINE

Instructions for Use

Description

The **SightLife Surgical Recipient Vacuum Trephine** has an ultra-sharp circular blade held inside an outer suction ring and attached to a spring-loaded syringe. Evacuating the ring when the trephine is placed in contact with the cornea, creates a vacuum which holds the trephine firmly in place on the cornea. Crosshairs aid centration over the optical axis.

Intended Use

The **SightLife Surgical Recipient Vacuum Trephine** is a single use trephine designed for full penetrating keratoplasty.

CAUTIONS

- The **SightLife Surgical Recipient Vacuum Trephine** blade is ultrasharp and penetrates the corneal tissue more efficiently than a traditional trephine. This trephine should not be used for partial lamellar keratoplasty.
- When 'setting' the blade at the beginning of the procedure, make certain that the blade is not in advance of the inner tube.
- Do not remove the blade completely from the body of the trephine as this may cause damage to the blade edge.
- Do not use if the sterile pack has been opened or damaged.

Surgical Technique

1. The trephine is supplied with one grooved marker pin on the capstan to aid blade control.
2. For safety reasons the trephine is supplied with the blade retracted. Before starting the procedure the blade needs to be set to the correct starting position. Using the microscope, set the blade by aligning the blade edge with the edge of the inner tube (Diagram 1).
3. When the blade is set, note the position of the grooved marker pin for future reference.

The 'Count back 8, Count forward 4' Principle

4. **Count back 8.** Retract the blade 8 quarter turns anticlockwise before positioning the trephine on the cornea. (Diagram 1)
5. Depress the plunger of the suction syringe and hold.
6. Centre the trephine correctly on the cornea using the crosshair.
7. Press down gently on the trephine to seat it on the cornea.
8. Release the plunger abruptly and wait approximately 10 seconds to allow the vacuum to become fully effective.
9. Lift the trephine gently to test whether the vacuum has taken and then allow the cornea to return to its normal curvature before rotating the blade .
10. **Count forward 4.** Advance the blade by 4 quarter turns to bring the blade to just above the surface of the cornea. The blade travels 63µm (0,063mm) for every quarter turn and 250µm (0,25mm) per complete rotation
11. Continue to rotate the blade slowly and gently clockwise to proceed with the cut.
12. Stop the blade rotation at the first sign of aqueous humour.
13. Release the suction by depressing the plunger of the syringe or by cutting the tubing and remove the trephine.
14. Complete the procedure in the desired manner.

Sterilisation

- The SightLife Surgical Recipient Vacuum Trephine is supplied sterile for SINGLE USE ONLY.
- DO NOT RE-STERILISE.

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:

1. Single use devices have not been validated for re-use.
If you re-use a device you may be held **Legally Liable** for the safe performance.
2. Cross-contamination and infection risks to patients . Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:
 - **Plastics:** Can be weakened, warped or become brittle.
 - **Metals:** Can be damaged or subjected to rusting. Blade edges will be damaged or blunt. Needles may be bent from the original shape.
 - **Adhesives:** Can degrade if re-sterilised causing parts to dislodge.
 - **Other materials:** May degrade, becoming unacceptable when compared to original manufacturing criteria.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

SightLife Surgical RECIPIENT VACUUM TREPHINE

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PRODUCT DETAILS

<i>SLS-10-650</i>	Recipient Vacuum Trephine 6.50 mm
<i>SLS-10-675</i>	Recipient Vacuum Trephine 6.75 mm
<i>SLS-10-700</i>	Recipient Vacuum Trephine 7.00 mm
<i>SLS-10-725</i>	Recipient Vacuum Trephine 7.25 mm
<i>SLS-10-750</i>	Recipient Vacuum Trephine 7.50 mm
<i>SLS-10-775</i>	Recipient Vacuum Trephine 7.75 mm
<i>SLS-10-800</i>	Recipient Vacuum Trephine 8.00 mm
<i>SLS-10-825</i>	Recipient Vacuum Trephine 8.25 mm
<i>SLS-10-850</i>	Recipient Vacuum Trephine 8.50 mm
<i>SLS-10-875</i>	Recipient Vacuum Trephine 8.75 mm
<i>SLS-10-900</i>	Recipient Vacuum Trephine 9.00 mm
<i>SLS-10-950</i>	Recipient Vacuum Trephine 9.50 mm

<i>Packaging</i>	1 per box Sterile SINGLE USE ONLY
------------------	---

ACCESSORIES

PRODUCT DETAILS

<i>Description</i>	Donor Trephine Punch with standard trephine Sizes 6.5mm - 9.50mm in 0.25mm increments
--------------------	--

<i>Packaging</i>	1 per box Sterile SINGLE USE ONLY
------------------	---

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