


# Reshaping Sight with CTAK

## Your Onboarding Guide

Below is an overview of our onboarding process and key responsibilities to ensure seamless integration into your practice.



### CTAK Onboarding Process

1. **Training:** Comprehensive wet lab sessions and online modules to ensure understanding of techniques.
2. **Patient Selection:** Guidance on identifying 3-4 suitable patients before scheduling a wet lab.
3. **Sending Scans:** Upload files to the secure CorneaGen Online Surgeon Portal or send via email to [CTAK@CorneaGen.com](mailto:CTAK@CorneaGen.com) for review.
4. **Surgical Plan:** Receive surgical plan options for each patient based on their specific needs.
5. **Order\*:** Place order using *Tissue Request Form* and include a copy of your selected surgical plan.
6. **Surgery\*\*:** Execute the procedure with support from CTAK's tools and resources.

*\*Please anticipate a 5-week period from CTAK order to delivery date.*

*\*\*To honor the processing of donor generated tissue and quality controls, any cancellation will result in invoiced responsibility to surgical facility.*

### Practice Responsibilities

To ensure a successful partnership, the following items are required from your practice:

- **Reimbursement:** Establish protocols for patient billing and insurance claims. CorneaGen will share our *CTAK Billing Primer* to support the process.
- **Tomography:**
  - **Pentacam:** 2012 or newer model with iris camera. HR, AXL, and AXL Wave.
  - **Galilei G4 & G6** with CTAK module (*software update released July 2025*).
    - **Files:** Securely upload files to the CorneaGen Online Surgeon Portal or send via email to [CTAK@CorneaGen.com](mailto:CTAK@CorneaGen.com)
- **Instruments:** CorneaGen partners with Corza for required surgical instrument for CTAK procedures. Your CorneaGen team will assist in aligning the local Corza representative with the designated surgical facility.
- **Laser:** Femto laser with corneal ring application installed – may require laser representative assistance. The following lasers are currently applicable for the CTAK procedure:
  - **Alcon** FS200 and LenSx, **Ziess** VisuMax 500, **Ziemer** Z4, Z6, Z8, **J&J** IntraLase (iFS)



Scan for additional CTAK reimbursement resources.



[Learn More](#)

To learn more visit [corneagen.com/CTAK](https://corneagen.com/CTAK) scan or call (800) 858-2020.

# Your Onboarding Guide

**OCULUS Pentacam General Overview**

**Corneal Topography Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Tomography Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Thickness Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Curvature Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Astigmatism Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Irregularity Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Elevation Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Thickness Map:**

The map shows the thickness of the cornea across the entire surface. The color scale ranges from blue (thin) to red (thick). The central thickness is approximately 550 micrometers.

**Corneal Curvature Map:**

The map shows the curvature of the cornea across the entire surface. The color scale ranges from blue (flat) to red (steep). The central curvature is approximately 43.80 D.

**Corneal Astigmatism Map:**

The map shows the astigmatism across the entire surface. The color scale ranges from blue (low) to red (high). The central astigmatism is approximately 0.00 D.

**Corneal Irregularity Map:**

The map shows the irregularity across the entire surface. The color scale ranges from blue (low) to red (high). The central irregularity is approximately 0.00 D.

**Corneal Elevation Map:**

The map shows the elevation across the entire surface. The color scale ranges from blue (low) to red (high). The central elevation is approximately 0.00 D.

**Corneal Thickness Profile:**

The profile shows the thickness of the cornea along a horizontal line. The thickness is approximately 550 micrometers.

**Corneal Curvature Profile:**

The profile shows the curvature of the cornea along a horizontal line. The curvature is approximately 43.80 D.

**Corneal Astigmatism Profile:**

The profile shows the astigmatism along a horizontal line. The astigmatism is approximately 0.00 D.

**Corneal Irregularity Profile:**

The profile shows the irregularity along a horizontal line. The irregularity is approximately 0.00 D.

**Corneal Elevation Profile:**

The profile shows the elevation along a horizontal line. The elevation is approximately 0.00 D.

**Corneal Thickness Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Curvature Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Astigmatism Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Irregularity Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Elevation Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Thickness Map:**

The map shows the thickness of the cornea across the entire surface. The color scale ranges from blue (thin) to red (thick). The central thickness is approximately 550 micrometers.

**Corneal Curvature Map:**

The map shows the curvature of the cornea across the entire surface. The color scale ranges from blue (flat) to red (steep). The central curvature is approximately 43.80 D.

**Corneal Astigmatism Map:**

The map shows the astigmatism across the entire surface. The color scale ranges from blue (low) to red (high). The central astigmatism is approximately 0.00 D.

**Corneal Irregularity Map:**

The map shows the irregularity across the entire surface. The color scale ranges from blue (low) to red (high). The central irregularity is approximately 0.00 D.

**Corneal Elevation Map:**

The map shows the elevation across the entire surface. The color scale ranges from blue (low) to red (high). The central elevation is approximately 0.00 D.

**Corneal Thickness Profile:**

The profile shows the thickness of the cornea along a horizontal line. The thickness is approximately 550 micrometers.

**Corneal Curvature Profile:**

The profile shows the curvature of the cornea along a horizontal line. The curvature is approximately 43.80 D.

**Corneal Astigmatism Profile:**

The profile shows the astigmatism along a horizontal line. The astigmatism is approximately 0.00 D.

**Corneal Irregularity Profile:**

The profile shows the irregularity along a horizontal line. The irregularity is approximately 0.00 D.

**Corneal Elevation Profile:**

The profile shows the elevation along a horizontal line. The elevation is approximately 0.00 D.

**Corneal Thickness Data:**


Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Curvature Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Irregularity	0.00 D

**Corneal Astigmatism Data:**

Parameter	Value
K1	43.80 D
K2	43.80 D
Kavg	43.80 D
Asigmatism	0.00 D
Ir	

[illegible][illegible]

CorneaGen<sup>TM</sup>