

Donor Trephine Punch

Precise. Accurate. Controlled.

SHARPEST BLADE

Ultra-thin, ultra-sharp blade to minimize undercut

HIGHEST VISIBILITY

High-visibility guide allows full visualization throughout the process to help accurately center and reposition the cornea

MOST CONTROL

Anatomically-correct bowl allows the cornea to center and stay in position, and finger grips allow for a stable base

Available with or without vacuum



For more information or to order,
please contact us at 844-526-7632 or
visit www.CorneaGen.com/DonorPunch

CorneaGen[™]

Donor Trephine Punch



Size Options

All sizes available with or without vacuum.

Trephine size (mm)							
6.00	6.25	6.50	6.75	7.00	7.25	7.50	7.75
8.00	8.25	8.50	8.75	9.00	9.25	9.50	9.75

Material Specification

Product Component	Specification
Trephine Blade	Stainless Steel
Trephine Body	Aluminum
Punch Base	Polycarbonate
Syringe with Luer	Polycarbonate, Luer-Lok (5ml capacity)
Outer Box/Carton	500 Micron Printed White Boxboard

Intended Use

The CorneaGen Donor Trephine Punch is designed specifically for preparing donor graft tissue. Each Trephine is designated for 'Single Use Only' and should not be resterilized by any method or reused.

Instructions for Use

The instructions for use and suggested surgical technique are supplied in the form of a multi-lingual leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are explained in each language. A leaflet is supplied with each product.

Sterilization

Products are sterilized by Ethylene Oxide (EO) according to a validated cycle and the requirements of BS EN ISO 11135-1:2007, by an approved sterilization sub-contractor, Synergy Health Plc. Doncaster, England. The products have a declared sterile shelf life of five years. Devices are designated for Single Use Only.

MANUFACTURED FOR CORNEAGEN