EndoSerter®-PL

FDA approved pre-loaded instrument for DSEK

DSEK REDEFINED

The exclusive FDA approved EndoSerter[®]-PL is a single-use graft insertion device pre-loaded with expertly processed Nano-Thin[™] or Ultra-Thin DSEK tissue. Successfully used in thousands of surgeries, the EndoSerter[®]-PL is the choice device for corneal surgeons.

EndoSeite

- \bigotimes
- Designed to allow consistent allograft positioning
- \bigotimes
- Simplifies the procedure, maximizing O.R. efficiency
- \bigcirc
- Reduces endothelium manipulation in the O.R.

EXPERIENCE DSEK REDEFINED

Contact us at (844) 526-7632 or visit www.CorneaGen.com



EndoSerter®-PL

Benefits

All the benefits of the EndoSerter[®] taken to the next level with pre-loaded DSEK tissue.

Simplifies DSEK Surgery

- Consistent allograft positioning
- Single 4mm incision
- Integral irrigation

Maximizes O.R. Efficiency

- Simplifies the procedure
- Saves valuable preparation time

Reduces Endothelial Cell Loss

- Minimizes environmental exposure
- Decreases endothelium manipulation in the O.R.

Specification Options

Customize your request with your preference for tissue thickness, age range, and cell count.

Tissue Thickness

- Nano-Thin[™] DSEK, less than 60µ
- Ultra-Thin DSEK, less than 100µ

Graft Size

≤ 8.0mm in allograft diameter

EndoSerter[®]-PL is validated for storage of donor tissue for up to a maximum of 48 hours.



101 N. Chestnut Street Ste 303 Winston-Salem, NC 27101 USA (844) 526-7632 CorneaGen.com SPS@CorneaGen.com

Request Your In-Service Today

EndoSerter- PL

An EndoSerter[®]-PL product in-service is required prior to initial use. Contact our team of expert Surgical Product Specialists to schedule your one-on-one in-service consultation.

Call (844) 526-7632 or email SPS@CorneaGen.com to get started.